



Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Product #:

2106028

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Powder reagent in sealed packages for the determination of Ortho-Phosphate. For 10mL samples; pk/1000.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/820
Method:	8048
Number of tests:	1000
Parameter:	Phosphate, Reactive
Platform:	Powder Pillows
Quantity:	1000 /pk
Range:	0.02 - 2.50 mg/L PO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)