



PhosVer 3 Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Product #:

212599

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagent in sealed package for the determination of Phosphorus, reactive (25mL sample; 100 tests).

Specifications

Instrument:	Best for use with DR2000, DR3000, DR900; or on other spectrophotometer with use of the 25-mL sample cell.
Method:	8048
Number of tests:	100
Parameter:	Phosphorus
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	0.02 - 2.50 PO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)