



Molybdovanadate Reagent Solution (100 mL MDB)

Product #:

2076032

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Molybdovanadate Reagent for the determination of ortho Phosphate, 0.3 - 45.0 mg/L PO₄

Specifications

Instrument:	All except DR/850, DR/820, PC II
Method:	Molybdovanadate
Number of tests:	50
Parameter:	Phosphate, ortho
Platform:	Solution
Quantity:	100 mL
Range:	0.3 - 45.0 mg/L PO ₄
Storage conditions:	10 °C to 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)