



## Low Range Orthophosphate Chemkey Reagents (box of 25)

|                     |                |
|---------------------|----------------|
| <b>Product #:</b>   | <b>8636600</b> |
| ZAR Price:          | Contact Hach   |
| Ships within 1 week |                |

Low Range Orthophosphate Chemkeys, range 0.20 - 4.00 mg/L PO<sub>4</sub>. For SL1000 Portable Parallel Analyser (PPA).

---

### Specifications

|            |                                  |
|------------|----------------------------------|
| Parameter: | Orthophosphate, LR               |
| Platform:  | Chemkeys                         |
| Quantity:  | pk/25                            |
| Range:     | 0.20 - 4.00 mg/L PO <sub>4</sub> |