



Phosphate, ortho, TNT, 0.06 - 5.00 mg/L PO₄

Product #:

2742545

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests.

Specifications

Method:	Ascorbic Acid
Number of tests:	25 - 50
Parameter:	Reactive (Orthophosphate)
Platform:	Test 'N Tube
Range:	0.06 - 5.00 mg/L PO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)