



Phosphate, total, TNT, 0.06 - 3.5 mg/L PO₄

Product #:

2742645

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Digestion Required:	Yes
Method:	8190
Number of tests:	25 - 50
Parameter:	Total Phosphorous
Platform:	Test 'N Tube
Range:	0.06 - 3.50 mg/L PO ₄
Reference method:	SM 4500 P-E
Special Remarks:	Set contains: Cuvette test tubes, Potassiumpersulfate Powder Pillows (2084766), PhosVer 3 Powder Pillows (2106046), Sodiumhydroxyde solution 1.54 N (2743042), Sulfuric acid solution 1.37 N (2742900)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)