



## Phosphate, total & ortho, TNT, 0.06 - 3.5 mg/L PO<sub>4</sub>

Product #:

2742745

ZAR Price:

Contact Hach

Ships within 3 days



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

---

### Specifications

Digestion Required:	Yes
Method:	Ascorbic Acid
Number of tests:	25 - 50
Parameter:	Total & Acid Hydrolyzable
Platform:	Test 'N Tube
Range:	0.06 - 5.00 mg/L PO <sub>4</sub>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)