



Phosphate, total, TNT, 1.0 - 100 mg/L PO_4

Product #:

2767245

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

| | |
|---------------------|---|
| Digestion Required: | Yes |
| Method: | 10127 |
| Number of tests: | 25 - 50 |
| Parameter: | Total Phosphorous |
| Platform: | Test 'N Tube |
| Range: | 1.0 - 100 mg/L PO_4 |
| Special Remarks: | Set contains: Cuvette test tubes, Potassiumpersulfate Powder Pillows (2084766), Molybdovanadate Reagent (2076026), Sodiumhydroxide solution 1.54 N (2743042), Sulfuric acid solution 1.37 N (2742900) |
| Storage conditions: | 10 - 25 °C |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)