



Ortho Phosphate, TNT, 1.0 - 100 mg/L PO₄

Product #:

2767345

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Method:	8114
Number of tests:	25 - 50
Parameter:	Reactive (Orthophosphate)
Platform:	Test 'N Tube
Range:	1.0 - 100.0 mg/L PO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)