



## COD Digestion Vials, Low Range, pk/150

**Product #:**

2125815

ZAR Price:  
Available

Contact Hach



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

**Contains mercury. Dispose of according to local, state and federal regulations.**

For determination of low range Chemical Oxygen Demand (COD) by the Reactor Digestion method.

US EPA approved for wastewater analysis using Hach Method 8000.

Range: 3-150 mg/L COD. Pack of 150 vials

### **EPA Approved**

Low Range (LR) and High Range (HR) vials are approved by the USEPA for NPDES reporting of wastewater COD concentrations.

### **Simpler and Safer**

Prepared COD reagents eliminate the need for measurement of corrosive and toxic chemicals. Just add sample to the COD vial, digest and measure. Photometric measurements provide accurate results without the equipment and effort required for titrations.

### **Easy Disposal**

Small reagent volumes used in Hach COD vials minimize the volume of spent silver, mercury and chromium waste. Customers in the contiguous 48 United States benefit from the availability of the ezCOD waste recycling program, which provides comprehensive transportation and disposal of COD vials.

---

## Specifications

|                     |                          |
|---------------------|--------------------------|
| Digestion Required: | Yes                      |
| EPA compliant:      | Yes                      |
| Instrument:         | All except DR/820, PC II |
| Method:             | 8000                     |
| Number of tests:    | 150                      |
| Parameter:          | Chemical Oxygen Demand   |
| Range:              | 3 - 150 mg/L COD         |

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)