



## Electrode Cleaning solution, RENOVO.N, for Clean Water Samples, 250 mL

**Product #:**

**S16M001**

**ZAR Price:**

Contact Hach

Ships within 1 week



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Electrode cleaning solution for pH sensing element regular cleaning maintenance. Convenient for soft cleaning of the pH sensing element after use in clean water samples. 250 mL bottle.

---

### Specifications

Description:	After measurements in aqueous solutions, cleaning with a mild solution such as RENOVO.N is sufficient.
Formula:	KATHON CG, DECONEX 11
Parameter:	Cleaning solution
Quantity:	250 mL
Solution Type:	Cleaning