



## GK ANNEX Maintenance Kit for Electrodes, Complete (Radiometer Analytical)

Product #:

S91M001

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

The GK ANNEX Maintenance Kit contains all you need for regular maintenance of combined pH, glass, and reference electrodes with saturated potassium chloride salt bridge. It contains: solutions for normal and intensive cleaning, KCl as saturated solution and crystals for refilling.

---

### Specifications

Description:	The Kit contains all you need for regular maintenance of combined pH, glass and reference electrodes with saturated potassium chloride salt bridge
Parameter:	GK ANNEX Electrode Maintenance Kit
Special Remarks:	Kit contains: Solutions for normal and intensive cleaning, KCl as saturated solution and crystals for refilling.