



## Electrode Cleaning Solution for Porous Pin/Diaphragm Junction, 250 mL

**Product #:**

**C20C380**

**ZAR Price:**

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Electrode cleaning solution for combined pH/Metal porous pin/diaphragm cleaning from Sulphides or AgCl precipitate/deposit/clogging.

---

### Specifications

Description:	Porous pin/diaphragm cleaning
Formula:	THIOUREA
Quantity:	250 mL
Solution Type:	Cleaning