



## Buffer Solution, pH 1.09, 500 mL

Product #:

S11M009

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis\*. Colourless pH buffer, pH is  $1.09 \pm 0.02$  at 25°C solution. 500 mL bottle

\*Available for download from [uk.hach.com/CoA](http://uk.hach.com/CoA): using product part number and lot number.

---

### Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
pH Value(s):	$1.09 \pm 0.02$ pH at 25°C
Quantity:	500 mL
Special Remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)