



Buffer Solution, pH 9.23, 500 mL

Product #: S11M011
ZAR Price: Contact Hach
Ships within 3-5 weeks

pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis*. Colourless pH buffer, pH is 9.23 ± 0.02 at 25°C solution. 500 mL bottle

*Available for download from uk.hach.com/CoA: using product part number and lot number.

Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
pH Value(s):	9.23 ± 0.02 pH at 25°C
Quantity:	500 mL
Special Remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)