



## Buffer Solution, pH 4.65, 500 mL

**Product #:** S11M010  
**ZAR Price:** Contact Hach  
**Available**

pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis\*. Colourless pH buffer, pH is  $4.65 \pm 0.02$  at 25°C solution. 500 mL bottle

\*Available for download from [uk.hach.com/CoA](http://uk.hach.com/CoA): using product part number and lot number.

---

### Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
pH Value(s):	$4.65 \pm 0.02$ pH at 25°C
Quantity:	500 mL
Special Remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)