



Crucible, Gooch, 40 mL Capacity

Product #:	2098641
ZAR Price:	Contact Hach
Ships within 1 week	

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Specifications

Capacity:	40 mL
Quantity:	1 /pk