



## Sartorius Biosart 100 Monitor

<b>Product #:</b>	<b>2521101</b>
ZAR Price:	Contact Hach
Call for ship date	

Sartorius Biosart 100 Monitors, Case of 48 Individually Packaged Monitors. The Sartorius Biosart 100 Monitors are incredibly unique and efficient. By drawing vacuum across the membrane, your sample is quickly and efficiently filtered through the included 0.45 µm filter. When you are done filtering, the monitor quickly converts into its own Petri dish, saving you time and money, while eliminating waste. Simply cap the bottom of the monitor, remove the cylinder side tube, add media, and cap the Petri dish. The Sartorius Biosart 100 Monitors are the most complete, quickest, and easiest to use filtering device that we have seen and are proud to offer this and other high-quality Sartorius branded items.

**Integrated Petri dish Innovative design offers sterile, disposable monitor that eliminates steps, saving time and money**

**Integrated membrane and absorbent pad**

**Integrated Petri dish**

**Also available in bulk**

---

### Specifications

Capacity:	100 mL
Material:	Cellulose Nitrate
Package Type:	Individual