



## Filter, TSS glass fiber, pore size 1.5 µm, diameter 30 mm, 100/pk

Product #: 2198730

 **Hazardous**

**OBsolete ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

---

### Specifications

Diameter:	30 mm
Material:	Glass Microfiber
Pore Size:	1.5 µm
Quantity:	100 /pk