



QL Filters, Pore Size 8 - 12 μ m, Diameter 42.5 mm, Circular, 100/pk, Qualitative

Product #: 50651

OBSOLETE ITEM
This item is no longer available.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifications

Diameter:	42.5 mm
Material:	Cellulose
Quantity:	100 /pk