



Filter, glass fiber, diameter 47 mm, preweighed standard methods TSS procedure, 100/pk

Product #: 2546100
ZAR Price: Contact Hach
Ships within 1 week

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Specifications

Diameter:	47 mm
Material:	Glass Microfiber
Pore Size:	1.5 μ m
Quantity:	100 /pk
Units:	100/pk