



Non-Sterile Nylon Syringe Filter, 33mm, 0.2um, 50/PK

Product #:	2513605
ZAR Price:	Contact Hach
Ships within 1 week	

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. Choose PVDF filters for solutions containing proteins. Choose Nylon filters for broad chemical compatibility and organic solutions. Standard Female Luer- Lok inlet with standard male slip Luer outlet. Hydrophilic.

Specifications

Diameter:	33 mm
Material:	Nylon
Pore Size:	0.20 μ m
Quantity:	50 /pk
Units:	50/pk