



Hach AE Grade Filters, 42.5 mm, 100/pk

Product #: 2551455

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hach Grade AE filters are comparable to Pall Type A/E filters and the 47 mm diameter is compatible with the magnetic filter holder #1352900. Filters feature fast flow rates, fine porosity and particle retention of 1.0µm which make them ideal for many dissolved and suspended solids testing procedures including Standard Methods 2540 C, 2540 D and 2540 E. This binder-free, borosilicate glass microfiber filter has a DOP efficiency of 99.98% (0.3µm) and a maximum operating temperature of 550 °C in air.

Specifications

Diameter:	42.5 mm
Material:	Glass Microfiber
Pore Size:	1.0 µm
Quantity:	100 /pk