



Filter, extra thick, 25 mm diameter, OriFlo filter holder, pk/100

Product #:

2676700

ZAR Price:
Available

Contact Hach

Designed for rapid, qualitative filtration of particulate matter from aqueous samples using a 25 mm glass fiber filter or other medium with porosity $\geq 1.0 \mu\text{m}$ (micron). Used in Hach Method 10014 for Ultra Low Range Chlorine in Wastewater (EPA Accepted method).

Specifications

Diameter:	25 mm
Quantity:	100 /pk