



## Hach AE Grade Filters, 150 mm, 100/pk

**Product #:** 2551463

**OBSOLETE ITEM**  
This item is no longer available.

Hach Grade AE filters are comparable to Pall Type A/E filters and the 47mm diameter is compatible with the magnetic filter holder #1352900. Filters feature fast flow rates, fine porosity and particle retention of 1.0µm which make them ideal for many dissolved and suspended solids testing procedures including Standard Methods 2540 C, 2540 D and 2540 E. This binder-free, borosilicate glass microfiber filter has a DOP efficiency of 99.98% (0.3µm) and a maximum operating temperature of 550 °C in air.

---

### Specifications

Diameter:	150 mm
Material:	Glass Microfiber
Pore Size:	1.0 µm
Quantity:	100 /pk