



Hach AE Grade Equivalent Filters, 47 mm, 100/pk

Product #:	2551456
ZAR Price:	Contact Hach
Available	

Hach Grade 162 filters are equivalent to Pall Type A/E filters. The 47mm diameter is compatible with magnetic filter holder PN 1352900. Filters feature fast flow rates, fine porosity and particle retention of 1.0 μm which make them ideal for many dissolved and suspended solids testing procedures including Standard Methods 2540 C, 2540 D and 2540 E. This binder-free, borosilicate glass microfiber filter has a DOP efficiency of 99.98% (0.3 μm) and a maximum operating temperature of 550 °C in air.

Specifications

Diameter:	47 mm
Material:	Glass Microfiber
Pore Size:	1.0 μm
Quantity:	100 /pk