



QL filters, pore size 8 - 12 μm , diameter 70 mm, circular, 100/ pk, qualitative

Product #: 50653

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered
hazardous under some shipping conditions.

If necessary, we will change your selected
shipping method to accomodate these items.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifications

Diameter:	70 mm
Material:	Cellulose
Quantity:	100 /pk