



## QL filters, pore size 8 - 12 $\mu$ m, diameter 125 mm, circular, 100/pk, qualitative

<b>Product #:</b>	<b>50657</b>
ZAR Price:	Contact Hach
Ships within 2 weeks	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Specifications

Diameter:	125 mm
Material:	Cellulose
Quantity:	100 /pk