



Filter, glass fiber, diameter 55 mm, preweighed standard methods TSS procedure, 100/pk

Product #:	2546120
ZAR Price:	Contact Hach
Available	

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Specifications

Diameter:	55 mm
Material:	Glass Microfiber
Pore Size:	1.5 μ m
Quantity:	100 /pk
Units:	100/pk