



Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Product #: 2546140
ZAR Price: Contact Hach
Ships within 3-5 weeks

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Specifications

Diameter:	90 mm
Material:	Glass Microfiber
Pore Size:	1.5 μm
Quantity:	100 /pk
Units:	100/pk