



## Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 400/pk

<b>Product #:</b>	<b>2546240</b>
ZAR Price:	Contact Hach
Call for ship date	

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

---

### Specifications

Diameter:	90 mm
Material:	Glass Microfiber
Pore Size:	1.5 $\mu$ m
Quantity:	400 /pk
Units:	400/pk