



Xenomax TSS-FVS Filter 42.5mm, 100/pk

Product #:	2511306
ZAR Price:	Contact Hach
Available	

Xenomax TSS-FVS Filters are specifically designed to ensure compliance with the filter preparation requirements of Standard Methods 2540 D and Standard Methods 2540 E (Fixed and Volatile Solids Ignited at 550°C). Independent 3rd party evaluation of Xenomax and 934-AH filters in both methods give demonstrably equivalent results.

In addition, these 1.5 µm binder-free, fast flowing, fine porosity borosilicate glass microfiber filters have been optimized for better wet strength/handleability and high temperature (550°C) blank stability. FVS Filters are available in popular diameter sizes and ready for use in a wide variety of labware including Buchner funnels and Gooch Crucibles.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

Specifications

Diameter:	42.5 mm
Material:	Glass Microfiber
Pore Size:	1.5 µm
Quantity:	100 /pk