



## Xenomax TSS-FVS Filter 125mm, 100/pk

<b>Product #:</b>	<b>2511312</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Xenomax™ TSS-FVS Raw Borosilicate Glass Fiber Filters are specifically designed to ensure compliance with the filter preparation requirements of Standard Methods 2540 D and Standard Methods 2540 E (Fixed and Volatile Solids Ignited at 550°C). Independent 3rd party evaluation of Xenomax and 934-AH® filters in both methods give demonstrably equivalent results.

In addition, these 1.5 µm binder-free, fast flowing, fine porosity borosilicate glass microfiber filters have been optimized for better wet strength/handleability and high temperature (550°C) blank stability. FVS Filters are available in popular diameter sizes and ready for use in a wide variety of labware including Buchner funnels and Gooch Crucibles.

*934-AH is a registered trademark of Whatman, Inc./GE Healthcare.*

---

### Specifications

Diameter:	125 mm
Material:	Glass Microfiber
Pore Size:	1.5 µm
Quantity:	100 /pk