



## Glass Fiber Filter, Diameter 70 mm, Preweighed for Standard Methods TSS procedure, 400/pk

Product #: 2546235

**OBSELETE ITEM**  
This item is no longer available.

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

---

### Specifications

Diameter:	70 mm
Material:	Glass Microfiber
Pore Size:	1.5 $\mu$ m
Quantity:	400 /pk
Units:	400/pk