



Filter Disc, Syringe, 1.2 μ m, Glass Fiber

Product #: 2599734

 **Hazardous**

OBSELETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Whatman 25mm disposable filters contain four filtration layers which help reduce blockage and increase volume throughput, allowing you to filter more of your sample in less time while increasing lab efficiency.

1.2 μ m GD/X. Non-sterile. Female Luer Lock inlet. 150 pk

Specifications

Diameter:	25 mm
Pore Size:	1.2 μ m
Quantity:	150 /pk