



## Syringe Filter, 1.5µm 25mm

<b>Product #:</b>	<b>2599738</b>
ZAR Price:	Contact Hach
Available	

Whatman 25mm disposable filters contain four filtration layers which help reduce blockage and increase volume throughput, allowing you to filter more of your sample in less time while increasing lab efficiency.

1.5µm 934-AH. Non-sterile. Female Luer Lock inlet. 150 pk

---

### Specifications

Diameter:	25 mm
Pore Size:	1.5 µm
Quantity:	150 /pk
What's included?:	Includes: 150 filters.

---

### What's included?

Includes: 150 filters.