



**Be Right™**



## Centrifuge Tube, Conical Bottom, 10/pk

<b>Product #:</b>	<b>2115839</b>
ZAR Price:	Contact Hach
Available	

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

---

### Specifications

Capacity:	15 mL
Length:	120 mm
Quantity:	10 /pk