



**Be Right™**



## **Flask, Erlenmeyer, glass w/screw cap, 6/pk, 500 mL**

<b>Product #:</b>	<b>2089679</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### **Specifications**

Capacity:	500 mL
Material:	Glass Screw Cap
Units:	6/pk