



## Flask, Erlenmeyer, glass w/screw cap, 6/pk, 1000 mL

Product #:

2089683

 **Hazardous**

### **OBSOLETE ITEM**

This item is no longer available.

Suggested replacements:

- 2089653

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### Specifications

Capacity:	1000 mL
Material:	Glass Screw Cap
Units:	6/pk