



Bottle, BOD, 300 mL, 6/pk

Product #:	62106
ZAR Price:	Contact Hach
Ships within 1 week	

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

Specifications

Body Material:	Glass
Capacity:	300 mL
Material:	Borosilicate Glass
Method:	8043
Package size UOM:	6pk
Quantity:	6 /pk