



## Flask, Erlenmeyer, polymethylpentene capacity 125 mL, 6/pk

<b>Product #:</b>	<b>2089873</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

---

### Specifications

Capacity:	125 mL
Material:	Polymethylpentene
Type:	Erlenmeyer
Units:	6/pk