



Flask, volumetric, polypropylene, 100 mL

Product #:	1406042
ZAR Price:	Contact Hach
Ships within 1 week	

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifications

Body material:	Polypropylene
Capacity:	100 mL
Certifications:	Not Applicable
Units:	1 each