



Flask, Volumetric, Polypropylene, 500 mL

Product #:	1406049
ZAR Price:	Contact Hach
Ships within 2 weeks	

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifications

Body material:	Polypropylene
Capacity:	500 mL
Certifications:	Not Applicable
Units:	1 each