



Flask, Volumetric, Polypropylene, 1000 mL

Product #:

1406053

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Polypropylene flask and screw cap. Durable with excellent chemical resistance. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifications

Body material:	Polypropylene
Capacity:	1000 mL
Certifications:	Not Applicable
Units:	1 each