



## Buret, automatic, class A, 5 mL

Product #: 2055037

**OBSOLETE ITEM**  
This item is no longer available.

Calibrated to ASTM E694 standard. NIST traceable, Class A buret with 70 mL reservoir. 5.0 mL buret capacity with 0.1 mL subdivision (tolerance  $\pm 0.01$  mL). Reservoir stopper and #2 PTFE stopcocks included.

---

### Specifications

|                   |                  |
|-------------------|------------------|
| Body material:    | Class A Glass    |
| Capacity:         | 5 mL             |
| Special Features: | Buret, Automatic |