



Cylinder, graduated, mixing, glass, 25 mL \pm 0.3 mL, 0.5 mL divisions, glass stopper #13

Product #: 189640
ZAR Price: Contact Hach
Ships within 1 week

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifications

Body material:	Glass
Capacity:	25 mL
Certifications:	Class B
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing