



Beaker, Berzelius glass, Tall Form, without spout, 200mL 12/pk

Product #:

2276175

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

200 mL Berzelius style (tall form) beakers from Kimax, are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Glass
Capacity:	200 mL
Material:	Glass, Berzelius, Tall Form
Units:	12 /pk