



Centifuge Tube PS, 15 mL, pack of 6

Product #:

2308306

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

These 15ml centrifuge tubes feature molded graduation markings, and are made from high quality polystyrene Conica bottoml. These tubes do not include caps. Gradations are 2.5, .5, 1.0, 2.5, 5, 10, 12, 15mL. Used in the LeadTrak Pocket Colorimeter Kit and the Ascorbic Acid Test Kit. 6/pk

Specifications

Capacity:	15 mL
Length:	130 mm
Quantity:	6 /pk