



Bottle, sampling, polycarbonate, 125 mL, 6/pk

Product #:	2324333
ZAR Price:	Contact Hach
Available	

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifications

Body Material:	Plastic
Capacity:	125 mL
Material:	PC
Package size UOM:	6pk
Quantity:	6 /pk