



Cylinder, Graduated, 1000 mL, Certified

Product #:	2636253
ZAR Price:	Contact Hach
Ships within 3-5 weeks	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Specifications

Body material:	Glass
Capacity:	1000 mL
Certifications:	Class A
Height:	465 mm
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric